
## Notice 598

### Echantillon de Germanium intrinsèque

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
Fig. 1

### Contents

1. Purpose
2. Description
3. Operation
4. Technical data
5. Experiment literature
6. Equipment list

#### 1. PURPOSE

The supporting plate 11807.00 is used to demonstrate the temperature dependence of the electrical conductivity of a pure (i.e. undoped) semiconductor. Unlike a metallic conductor or a heavily doped

semiconductor, in which extrinsic conduction predominates, an undoped semiconductor exhibits a sharp rise in conductivity with rising temperature due to the increasing concentration of thermally produced electrons and holes (intrinsic conduction).

#### 2. DESCRIPTION

Fig. 1 shows the supporting plate with the individual functional units:

Semiconductor crystal 1 is an undoped germanium crystal with the dimensions 20 x 10 x 1 mm.

A suitable DC voltage is applied to sockets 2 through a series resistance which limits the current to a maximum of 30 mA (see Operation). These sockets are also used when measuring the voltage drop through the crystal.

The two plugs 3 on the back of the plate not only give mechanical support to the plate, e.g. in the multisocket distributor 06024.00, but also act as leads for the heating current (heating time limited, see Operation).

The copper-constantan thermocouple 4 provides a thermoelectric voltage of approx. 40  $\mu\text{V/K}$  which can be taken off through the two sockets 5.

### 3. OPERATION

#### Caution:

Germanium crystals are very brittle and consequently fragile. To prevent any bending of the supporting plate during assembly, the following points should be observed:

To affix the plate to the multisocket distributor 06024.00, the two plugs 3 are pushed into the distributor sockets by pressing on the front side of the plate between the nuts securing these plugs. To remove the plate, the back of the plate near the plugs is pressed with two fingers and pulled away from the distributor.

When wiring up the individual sockets, these are each given support by pressing against the back of the plate with the hand.

It is further recommended that the original packing should be used when storing the supporting plate.

In the experiment set-up, the multisocket distributor is mounted on a suitable stand. The two plugs 3 can also be used to take the heating current of 5 A at 6 V DC or AC via the distributor; however, the voltage source should not be switched on until just before the measurement is made.


#### Caution:

To prevent overheating, the heating current is switched off as soon as the thermoelectric voltage (sockets 5) reaches a value of 5 mV.

An appropriate DC voltage (maximum 42 V) is applied to sockets 2 via a protective resistor which prevents the control current from rising above 30 mA even when, at high temperatures, the resistance of the crystal is very low. To determine the electrical resistance of the crystal, the voltage drop in the crystal is measured across these sockets. In addition, the current flowing through the crystal is measured with a second meter.

To measure the thermoelectric voltage, a DC voltmeter (range multiplier 10 mV, or if need be 30 mV) is applied across sockets 5.

Fig. 2

The voltage drop in the crystal can be recorded as a function of the thermoelectric voltage by using an xy recorder in place of the two voltmeters.

Fig. 2 shows the principle of the measuring circuit used. The resistance was designed for a working voltage of approx. 18 V in order to get a current of just 1 mA. The current obtained is largely unaffected by the temperature because of the relatively large series resistance.

#### 4. TECHNICAL DATA

##### Crystal:

Dimensions 20 x 10 x 1 mm  
Material germanium, undoped  
Resistivity  $> 15 \Omega \cdot \text{cm}$   
Max. permissible current 30 mA

##### Other data:

Maximum temperature  $175^{\circ}\text{C}$

Thermocouple Cu-CuNi

Thermoelectric voltage coefficient approx.  $40 \mu\text{V/K}$

Heating current 5 A at 6 V AC or DC

#### 5. EXPERIMENT LITERATURE

PHYWE Experimental Units in Physics,  
Solid-state physics (in preparation)

#### 6. EQUIPMENT LIST

The following list gives the most important equipment used in connection with the supporting plate "Intrinsic conduction of germanium" 11807.00.

Order No.	Description
-----------	-------------

For power supply:

11704.93	Power supply, universal
or	
11706.93	Power supply, stabilized (heating current 5 A at 6 V AC or DC then required as an addition)

For measurement of both voltages:


11408.97	myt Recorder
or	
suitable voltmeters (3 V DC, 1 V DC, 300 mV DC and 10 mV DC) and, specially for measuring thermoelectric voltages:	

11100.00	Moving-coil instrument with
11110.01*	Range multiplier, amplifier 30 $\mu\text{V}$ to 10 mV
or	
11761.93	Measuring amplifier (with suitable display instrument)

For holding the supporting plate:

06024.00	Multisocket distributor
----------	-------------------------

\* Four 9 V batteries 07496.10 are required for this range multiplier

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
## Operating instructions

### 1 PURPOSE AND CHARACTERISTICS

The carrier plate allows to demonstrate the temperature dependence of the electric conductivity of a pure, that is, non doped, semiconductor. Contrarily to a metallic conductor, or to a strongly doped semiconductor, in which extrinsic conduction is dominant, in the case of a non doped semiconductor, a strong increase of conductivity is observed as temperature increases, due to an increasing concentration of thermal electrons and holes (intrinsic conductivity).

### 2 DESCRIPTION AND HANDLING

#### 2.1 Functional and control elements

A continuous current source is connected to the semiconductor crystal 1 over the sockets 2, whose current intensity must be limited by a protective resistor inserted before the crystal. The sockets 2 are also used to connect a voltmeter, to determine the drop of voltage across the crystal. The copper-constantan thermocouple 4 provides a thermoelectric voltage which can be tapped at the pair of sockets 5.

Together with the corresponding PHYWE distributor (cf. list of accessories), the pair of plugs 3 at the back side of the plate is used both for mechanical support of the carrier plate and to supply heating current.

#### 2.2 General set-up

Attention! Germanium crystals are very brittle and thus can easily break. The following recommendations should be followed in order to avoid bending of the carrier plate and thus mechanical stress to the crystal:

- when plugging the carrier plate into the corresponding distributor, the pair of plugs 4 is pushed into the distributor sockets, pressing onto the front side of the plate, between the fastening nuts. To remove the plate, a finger is inserted behind the plate near the plug and the latter is carefully pulled off the distributor.
- When connecting the single sockets, they are supported with the hand at the back of the plate.

A continuous voltage suited to the protective resistor is applied to the sockets 2 over a protective resistor, which is required to avoid an increase of the current intensity above 30 mA (especially if the crystal resistance drops to very low values at high temperatures). To heat the crystal, a maximum heating voltage of 5 V / 5 A is applied to the two sockets 5 or over the distributor. The heating voltage must be switched on only immediately before performing the measurement.

Temperature control is carried out using the integrated thermocouple and connecting a sensitive voltmeter (measuring range 10-30 mV).

**Attention!** As soon as the thermoelectric voltage has reached 5 V (heating time about 2 minutes), the heating current must be switched off, to avoid overheating of the carrier plate. A possible brown discoloration of the carrier plate, which may occur after heating to the maximum admissible temperature limit, has no effect on the function of the unit.

#### 2.3 Hints for the performance of the experiment

The distributor must be held with adequate support material for the experimental set-up.

As an alternative, an xy plotter may be used to record the drop of voltage over the crystal as a function of the thermoelectric voltage.

The principle of the used measurement circuit is represented in fig. 2. The 22 k $\Omega$  protective resistor is designed to stand an operating voltage of about 18 V, the used sample current being kept just below 1 mA on purpose. The relatively large protective resistor allows to keep the current intensity practically independent of temperature.


Fig. 2

### 3 TECHNICAL SPECIFICATIONS

Crystal material	germanium, non doped
Crystal dimensions	(20 x 10 x 1) mm
Specific resistivity	> 15 $\Omega$ cm
Maximum current	30 mA
Thermocouple	Cu-CuNi
Thermoelectric voltage coefficient	approx. 40 mV/K
Heating voltage	6 V $\equiv$ / approx. 5 A

### 4 LITERATURE REFERENCES

University Laboratory Course in Physics 1-3 16502.02

### 5 LIST OF ACCESSORIES

Distributor 06024.00

Universal power supply 13500.93

or

Controlled power supply 0-12 V/2 A 13505.93

Furthermore, an ammeter and a voltmeter are required.